



Parts list			
Qty	Part #	Description	Size
4	AF24L	Base Cabinet, X= No Toe Kick, (4 ea) 1/8 Drawers at Top & (1 ea) 1/4 Drawer at Bottom, (APP)	24 1/8"L x 29 3/4"H x 22"D
4	B47L	Wall Cabinet, X= 24" Height	47" L x 30" H x 13" D
2	CTFR-100	Flex Bench Counter Frame, 100" Long	100" x 23"
4	LBRACE	Flex Bench Lower Cross Brace	"As Shown"
6	LCTB	Flex Bench Counter Support Bracket	8" x 15"
6	POSTLEG	Flex Bench Side Leg w/ Post Standard	84" Tall
4	UHR-47	Upper Cross Brace, 47" Long	47" Long
8	BEBWC	Book End Bracket For Wall Cabinet Installation	"As Shown"
4	PM4812	WireMold Two Piece Steel Raceway	4' Long w/ Plugs at 12" O.C.
2	TRESPA	3/4" Thick Athlon Grade Trespa Countertop	"As Shown"

HANSON LAB FURNITURE, INC.
 814 MITCHELL RD.
 NEWBURY PARK, CA 91320
 PH: 805/498-3121
 FAX: 805/498-1855



HEWLETT-PACKARD
 PALO ALTO, CA.

DRAWN BY	DATE	SCALE	QUOTE NO.	REV. NO.
JDP	11/01/04	3/8" = 1'	4124	
SHEET NO:	JOB NO:			
1 OF 2				