



| TABLE # | LENGTH | WIDTH |
|-----------|--------|-----------|
| 30" X 36" | 34" | 22" & 28" |
| 30" X 48" | 46" | 22" & 28" |
| 30" X 54" | 52" | 22" & 28" |
| 30" X 60" | 58" | 22" & 28" |
| 30" X 72" | 70" | 22" & 28" |

- NOTES:
1. STANDARD FINISH IS EPOXY POWDER COATED STEEL
 2. ALSO AVAILABLE IN STAINLESS STEEL
 3. CASTER KITS AVAILABLE FOR MOBILE TABLES
 4. OVERHEAD SHELF AVAILABLE

ALL TOLERANCES ARE +/- .030 UNLESS OTHERWISE NOTED

| | | |
|---------------------|-----------------|-----------------------|
| DRAWN BY JDP | DATE 1/12/07 | SCALE 1/2"=1'-0" |
| SHEET NO: 1 OF 1 | JOB NO: | QUOTE NO: REV. NO. |

ADJ. HT. TABLE



HANSON LAB FURNITURE, INC.
814 MITCHELL RD.
NEWBURY PARK, CA 91320
PH: 805/498-3121
FAX: 805/498-1855