

The Science  
of  
Laboratory  
Construction



**Hanson**

## How do you build a better lab? Ask a scientist.

*"We don't just fill space; we build functional workstations around your exact process."*  
 Joe Matta, VP, Business Development

Ask the science professionals who founded and run Hanson Lab Furniture. Ask the Hanson lab planners – knowledgeable professionals who have the science and construction backgrounds to

- listen to what you say,
- understand the challenges you face,
- design according to your needs,
- deliver a complete lab, from plumbing to countertops, that works with your exact process – even if it means customizing or designing furniture just for you.

Most importantly, ask the thousands of scientists who have relied on Hanson to design, build, and deliver complete, high-quality labs for over 30 years.

The Hanson difference starts with the Hanson Lab Planner – a science professional who will stay with your



*Finishes shown: Wild Rose and Sandstone. Custom finishes available.*

project from your first meeting to the day you begin work in your new lab.

Your Lab Planner will recommend composite materials compatible with your chemistries, design and build functional workstations around your exact processes, and work with your Hanson Project Manager to coordinate every aspect of your project and to keep you informed every step of the way.



*Finish shown: Antique White with Silver Grey Trespa Countertops.*



*Standard Fume hoods are available in sizes from 3' to 10'. Custom sizes also available.*

## The secret to outstanding quality, high value, and the fastest installation times in the industry? Vertical integration.

*"As the Manufacturer, we can deliver your lab faster and to your exact specifications."*  
 Mike Hanson, President & CEO

We build all our own lab furniture, from casework to countertops. Complete control over the manufacturing process, and direct sales to our customers, mean that we can offer the highest quality at a competitive price. It also gives you the flexibility to build or customize solutions to your specifications.

Our steel casework, vented chemical fume hoods, and specialty safety cabinets are built by hand, by artisans with an average of more than ten years experience in metal fabrication. We use high-end stainless steel and first-rate



composite counters. Our quality details raise our products to a level of appearance, functionality, and durability you won't find anywhere else.

And because we're also general contractors, we can interface with facilities managers, subcontractors and engineers to keep your project on track while you concentrate on what you do best.




*Finishes shown: Sky Blue with Black Satin Epoxy countertops.*



# Quality. Durability. Functionality. And great looks, too.

You'll be spending a lot of time in your new lab – shouldn't it be great-looking, too?

We think so. That's why Hanson cabinets and countertops are attractively styled, and available in a wide variety of colors and finishes. (Our finishes are so beautiful, in fact, that some owners of custom homes have demanded our lab furniture to equip their gourmet kitchens.)



Yet they're durable enough to stand up to the rigors of real lab work and harsh chemicals, and to look great for years to come. Your Lab Planner can help you choose the finishes that are best for your chemistries.

Whether you work in research & development, pharmaceuticals, biomedical, or testing, the result is a Hanson lab that's functional, comfortable, safe, and science-smart.

## At Hanson, details make the difference.

### Hanson



**Hanson:** Vertically-welded supports provide greater load-bearing capacity and assure proper door and drawer alignment. They are ground flat and powder-coated, for corrosion resistance and long, attractive life.

### Others



**Others:** Non-welded components are more susceptible to deterioration and can go out of adjustment over time.



*Casework, shelving, countertops and vent hoods  
all designed and manufactured by Hanson.*

## The Hanson advantage.

- Direct-from-manufacturer value
- Competitively priced
- Ability to customize standard casework and fume hoods
- Complete line of lab furniture and accessories:
  - vented hoods
  - cabinets
  - sinks
  - faucets
  - gas valves
  - countertops
  - safety equipment
- 100% epoxy powder coat finish
- UL listed fume hoods
- Custom colors
- Complete refurbishing
- The fastest lead times in the industry



**Hanson**

The Science of  
Laboratory Construction

Hanson Lab Furniture, Inc.

814 Mitchell Road

Newbury Park, CA 91320

Contractor License # B799172

Phone: 805-498-3121

Fax: 805-498-1855

[www.hansonlab.com](http://www.hansonlab.com)